



Par Scientific

Apollo version 1.0.1.14

Dear user.

This is the version for Apollo 5 / Apollo 10 for reading the values from the dosimeter.

Viruscheck

The computer generating the CD has been continuously checked for Virus, and a scan has been made using Norton LifeLock 22.23.6.5 of 2023-09-05

Requirements:

.Net 4.72 or newer

Installation:

The software will only replace program files and leave all setup files alone. All settings and data stored will be kept unchanged.

If questions should arise please feel free to contact your local distributor or PAR Scientific directly.

PAR Scientific
Falstervej 8D
5800 Nyborg
Denmark

Phone: +45 63 12 90 00
E-mail: tech@ParScientific.com
<http://www.ParScientific.com>



Par Scientific

Last update: 2023-09-05

1.0.1.14 Added curve and histogram views
1.0.0.0 Porting program from DOS